

Data Sheet

HUMAN CHEMOKINE CC-2 (aa 78-92)

ANTIBODY, POLYCLONAL

| | |
|----------------------------|--|
| Catalog no.: | AE1115.1 / AE1115.2 |
| Immunogen: | Synthetic human Chemokine CC-2 (aa 78-92) poly Lysine conjugated (PGVQDCMKLKPYSI) |
| Synonyms: | C-C motif chemokine 15, Small-inducible cytokine A15 |
| Swiss-Prot No: | Q16663 |
| Gene Information: | Gene Name: CCL15 GeneID: 6359 |
| Host: | Rabbit |
| Matrix: | Serum |
| Specificity: | Human Chemokine CC-2 (aa 78-92) |
| Contents: | 20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest. |
| Known applications: | ELISA (1:1000) This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls. |
| Store at: | 2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided |
| Last updated on: | 30 June 2016 |

For research use only

Publishing research using AE1115? Please let us know so that we can cite your publication as a reference.

