

Data Sheet

HUMAN CHEMOKINE CC-1/3 (aa 21-35)

ANTIBODY, POLYCLONAL

Catalog no.:	AE 1117.1 / AE 1117.2
Immunogen:	Synthetic human Chemokine CC-1/3 (aa 21-35) poly Lysin- conjugated (TTYKIPRQRIMDYEE)
Synonyms:	C-C motif chemokine 14, Small-inducible cytokine A14
Swiss-Prot No:	Q16627
Gene Information:	Gene Name: CCL14 Gene ID: 6358
Host:	Rabbit
Matrix:	Serum
Specificity:	Synthetic human Chemokine CC-1/3 (aa 21-35)
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	ELISA (1:1,000) This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	19 August 2016

For research use only

Publishing research using AE1117? Please let us know so that we can cite your publication as a reference.

