



Data Sheet

HUMAN CHEMOKINE CC-1/3 (aa 21-35)

SYNTHETIC

Catalog no.:	AE1117AG.1 / AE1117AG.2
Name:	Human Chemokine CC-1/3 (aa 21-35)
Synonyms:	C-C motif chemokine 14, HCC-1/HCC-3, HCC-1(1-74), NCC-2, Small-inducible cytokine A14
Sequence:	TTYKIPRQRIMDYEE
MW:	1.98 kDa
Swiss-Prot No:	Q16627
Gene Information:	Gene Name: CCL14, NCC2, SCYA14 Gene ID: 6358
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration.
Store at:	2-8 °C (lyophilized); -20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	23 April 2018

For research use only

Publishing research using AE1117AG? Please let us know so that we can cite your publication as a reference.

