



DATA SHEET

HUMAN CHEMOKINE CC-1/3 (aa 21-35)

Catalog no.:	AE1117AG.1 / AE1117AG.2
Protein:	Human Chemokine CC-1/3 (aa 21-35) (TTYKIPRQRIMDYEE)
Reference no.:	LocusID 6358
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg/ 1 mg (lyophilized) Resuspend in 0,1% acetic acid for required concentration
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided

Only for research use

N:\WINWORD\DAT\DATEN\LAETTER\PROTEINE\A 1005AG-proBNP1-76rec.DOC