

Data Sheet

HUMAN CHEMOKINE CC-1/3 (aa 8-15)

ANTIBODY, POLYCLONAL

Catalog no.:	AE 1118.1 / AE 1118.2
Immunogen:	Synthetic human Chemokine CC-1/3 (aa 8-15) poly Lysin- conjugated (RGPYHPSE)
Synonyms:	C-C motif chemokine 14, Small-inducible cytokine A14
Swiss-Prot No:	Q16627
Gene Information:	Gene Name: CCL14 GeneID: 6358
Host:	Rabbit
Matrix:	Serum
Specificity:	Synthetic human Chemokine CC-1/3 (aa 8-15)
Contents:	20 µl / 100 µl (lyophilized) resuspend in 20 µl / 100 µl aqua bidest
Known applications:	ELISA (1:1,000) <small>This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.</small>
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided

For research use only

