



Data Sheet

HUMAN CHEMOKINE CC-1

ANTIBODY, POLYCLONAL

Catalog no.:	AE 1119.1 / AE 1119.2
Immunogen:	Native human Chemokine CC-1, poly Lysin conjugated
Synonyms:	None
Swiss-Prot No:	Q76BR7
Gene Information:	Gene Name: CCL1 GenelD: 6346
Host:	Chicken
Matrix:	IgY
Specificity:	Human Chemokine CC-1
Contents:	10 µg / 100 µg (lyophilized) resuspend in 10 µl / 100 µl aqua bidest
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided

For research use only

