

Data Sheet

HUMAN CHEMOKINE CC-1

ANTIBODY, POLYCLONAL

Catalog no.:	AE 1119.1 / AE 1119.2
Immunogen:	Native human Chemokine CC-1, poly Lysin conjugated
Synonyms:	None
Swiss-Prot No:	Q76BR7
Gene Information:	Gene Name: CCL1 GeneID: 6346
Host:	Chicken
Matrix:	IgY
Specificity:	Human Chemokine CC-1
Contents:	10 µg / 100 µg (lyophilized) Resuspend in 10 µl / 100 µl aqua bidest.
Known applications:	ELISA <small>This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.</small>
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	12 January 2018

For research use only

Publishing research using AE1119? Please let us know so that we can cite your publication as a reference.

