



DATA SHEET

HUMAN CHEMOKINE CC-1, NATIVE

Catalog no.:	AE1119AG.1 / AE1119AG.2
Protein:	Human Chemokine CC-1, native
Reference no.:	LocusID 6346
Source:	
Purity:	> 95% (HPLC)
Contents:	100 µg/ 1 mg (lyophilized) Resuspend in aqua bidest for required concentration
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided

Only for research use

N:\WINWORD\DAT\DATEN\LAETTER\PROTEINE\A 1005AG-proBNP1-76rec.DOC