



Data Sheet

HUMAN CHEMOKINE CC-1 SYNTHETIC

Catalog no.:	AE1119AG.1 / AE1119AG.2
Name:	Human Chemokine CC-1
Synonyms:	C-C motif chemokine 1, Small-inducible cytokine A1, T lymphocyte-secreted protein I-309
MW:	8.49 kDa
Swiss-Prot No:	P22362
Gene Information:	Gene Name: CCL1, SCYA1 Gene ID: 6346
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in aqua bidest. for required concentration.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	20 April 2018

For research use only

Publishing research using AE1119AG? Please let us know so that we can cite your publication as a reference.

