



Data Sheet

HUMAN CHEMOKINE CC-2

ANTIBODY, MONOCLONAL

Catalog no.:	AE 1120.1 / AE 1120.2
Immunogen:	Native human Chemokine CC-2, poly Lysine conjugated
Synonyms:	C-C motif chemokine 15, Small-inducible cytokine A15
Swiss-Prot No:	Q16663
Gene Information:	Gene Name:CCL15 GeneID: 6359
Host:	Chicken
Matrix:	IgY
Specificity:	Native human Chemokine CC-2
Contents:	10 µg / 100 µg (lyophilized) resuspend in 10 µl / 100 µl aqua bidest
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided

For research use only

