

Data Sheet

HUMAN CHEMOKINE CC-2

ANTIBODY, POLYCLONAL

Catalog no.:	AE 1120.1 / AE 1120.2
Immunogen:	Native human Chemokine CC-2, poly Lysine conjugated
Synonyms:	C-C motif chemokine 15, Small-inducible cytokine A15
Swiss-Prot No:	Q16663
Gene Information:	Gene Name:CCL15 Gene ID: 6359
Host:	Chicken
Matrix:	IgY
Specificity:	Native human Chemokine CC-2
Contents:	10 µg / 100 µg (lyophilized) Resuspend in 10 µl / 100 µl aqua bidest.
Known applications:	ELISA This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	12 January 2018

For research use only

Publishing research using AE1120? Please let us know so that we can cite your publication as a reference.

