



## *DATA SHEET*

# HUMAN CHEMOKINE CC-2, NATIVE

<b>Catalog no.:</b>	AE1120AG.1 / AE1120AG.2
<b>Protein:</b>	Human Chemokine CC-2, native
<b>Reference no.:</b>	LocusID 6359
<b>Source:</b>	
<b>Purity:</b>	> 95% (HPLC)
<b>Contents:</b>	100 µg/ 1 mg (lyophilized) Resuspend in 0,1% acetic acid for required concentration
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided

**Only for research use**

N:\WINWORD\DAT\DATEN\LAETTER\PROTEINE\A 1005AG-proBNP1-76rec.DOC