

DATA SHEET

HUMAN CHOLECYSTOKININ-RECEPTOR [CCK-R] TYPE B (aa 191-204)

Catalog no.:	AE 1127AG.1 / AE 1127AG.2
Immunogen:	Human CCK-R Type B (aa 191-204) (VYTVVQPVGPRVLQ)
Reference no.:	LocusID 887
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in 0,1% acetic acid for required concentration
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
References:	1. Kulaksiz H, Arnold R, Goke B, Maronde E, Meyer M, Fahrenholz F et al. Expression and cell-specific localization of the cholecystokinin B/gastrin receptor in the human stomach. Cell Tissue Res 2000;299:289-98

Only for research use

N:\WINWORD\DAT\DATEN\LAETTER\ANTIKOERPER\MISCELLANEOUS\AE 1127.DOC