

## *DATA SHEET*

# **HUMAN CHOLECYSTOKININ-RECEPTOR [CCK-R] TYPE B (aa 206-219)**

<b>Catalog no.:</b>	AE 1128AG.1 / AE 1128AG.2
<b>Immunogen:</b>	Human CCK-R Type B (aa 206-219) (VHRWPSARVRQTWS)
<b>Reference no.:</b>	LocusID 887
<b>Source:</b>	Synthetic
<b>Purity:</b>	> 95% (HPLC)
<b>Contents:</b>	100 µg / 1 mg (lyophilized) Resuspend in 0,1% acetic acid for required concentration
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
<b>References:</b>	1. Kulaksiz H, Arnold R, Goke B, Maronde E, Meyer M, Fahrenholz F et al. Expression and cell-specific localization of the cholecystokinin B/gastrin receptor in the human stomach. Cell Tissue Res 2000;299:289-98

**Only for research use**

N:\WINWORD\DAT\DATENBLAETTER\ANTIKOERPER\MISCELLANEOUS\AE 1128.DOC