



Data Sheet

HUMAN CHOLECYSTOKININ-RECEPTOR TYPE A [CCK-AR] (aa 180-195)

ANTIBODY, POLYCLONAL

Catalog no.:	AE1129.1 / AE1129.2
Immunogen:	Synthetic human CCK-R Type A (aa 180-195) polyLysine conjugated (SNLVPFTKNNNQTANM)
Synonyms:	Cholecystokinin-1 receptor (CCK1-R)
Swiss-Prot No:	P32238
Gene Information:	Gene Name: CCKAR, CCKRA GenelD: 886
Host:	Rabbit
Matrix:	Serum
Specificity:	Human CCK-R Type A (aa 180-195)
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	20 December 2011

For research use only

