



Data Sheet

HUMAN CHOLECYSTOKININ-RECEPTOR TYPE A [CCK-R] (aa 180-195)

SYNTHETIC

| | |
|--------------------------|---|
| Catalog no.: | AE1129AG.1 / AE1129AG.2 |
| Name: | Human CCK-R Type A (aa 180-195) |
| Synonyms: | CCK-A receptor (CCK-AR), Cholecystokinin-1 receptor (CCK1-R) |
| Sequence: | SNLVPFTKNNNQTANM |
| MW: | 1.79 kDa |
| Swiss-Prot No: | P32238 |
| Gene Information: | Gene Name: CCKAR, CCKRA Gene ID: 886 |
| Source: | Synthetic |
| Purity: | > 95% (HPLC) |
| Contents: | 100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration. |
| Store at: | 2-8 °C (lyophilized); -20 °C (dissolved) Repeated thawing and freezing must be avoided |
| Last updated on: | 23 April 2018 |

For research use only

Publishing research using AE1129AG? Please let us know so that we can cite your publication as a reference.

