



DATA SHEET

ANTI HUMAN CHOLECYSTOKININ-RECEPTOR [CCK-R] TYPE A (aa 197-210) ANTIBODY, POLYCLONAL

Catalog no.:	AE 1130.1 / AE 1130.2
Immunogen:	Synthetic human CCK-R Type A (aa 197-210) poly Lysin conjugated (RFLPPNDVMQQSWH)
Reference no.:	LocusID 886
Host:	Rabbit
Matrix:	Serum
Specificity:	Human CCK-R Type A (aa 197-210)
Contents:	20 µl / 100 µl (lyophilized) resuspend in 20 µl / 100 µl aqua bidest
Known applications:	Westernblot (1:2000) This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided

Only for research use

N:\WINWORD.DAT\IMMUNOCHEMIKALIEN\DATEN\BLAETTER\ANTIKOERPER\MISCELLANEOUS\AE1130.DOC