



Data Sheet

HUMAN CHOLECYSTOKININ-RECEPTOR [CCK-R] TYPE A (aa 197-210)

SYNTHETIC

Catalog no.:	AE1130AG.1 / AE1130AG.2
Name:	Human CCK-R Type A (aa 197-210)
Synonyms:	CCK-A receptor (CCK-AR), Cholecystokinin-1 receptor (CCK1-R)
Sequence:	RFLPPNDVMQQSWH
MW:	1.87 kDa
Swiss-Prot No:	P32238
Gene Information:	Gene Name: CCKAR, CCKRA Gene ID: 886
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration.
Store at:	2-8 °C (lyophilized); -20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	23 April 2018

For research use only

Publishing research using AE1130AG? Please let us know so that we can cite your publication as a reference.

