

## *DATA SHEET*

# HUMAN CHOLECYSTOKININ-RECEPTOR [CCK-R] TYPE A (aa 197-210)

<b>Catalog no.:</b>	AE 1130AG.1 / AE 1130AG.2
<b>Immunogen:</b>	Human CCK-R Type A (aa 197-210) (RFLPPNDVMQOSWH)
<b>Reference no.:</b>	LocusID 886
<b>Source:</b>	Synthetic
<b>Purity:</b>	> 95% (HPLC)
<b>Contents:</b>	100 µg / 1 mg (lyophilized) Resuspend in 0,1% acetic acid for required concentration
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided

**Only for research use**

N:\WINWORD\DAT\DATENBLAETTER\ANTIKORPER\MISCELLANEOUS\A 6100 2.DOC