

Data Sheet

HUMAN CYCLIC AMP [cAMP]

ANTIBODY, POLYCLONAL

Catalog no.:	AE1131.1 / AE1131.2
Immunogen:	Synthetic human succinylated cAMP, BSA conjugated
Host:	Rabbit
Matrix:	Serum
Specificity:	Synthetic human succinylated cAMP
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	RIA (1: 1000) This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	7 February 2018

For research use only

Publishing research using AE1131? Please let us know so that we can cite your publication as a reference.

