



DATA SHEET

ANTI HUMAN PNEUMADIN ANTIBODY, POLYCLONAL

Catalog no.:	AE1134.1 / AE1134.2
Immunogen:	Synthetic human Pneumadin (AGEPKLDAGV), poly Lysin conjugated
Reference no.:	P22103
Host:	Rabbit
Matrix:	Serum
Specificity:	Synthetic human Pneumadin
Contents:	20 µl / 100 µl (lyophilized) resuspend in 20 µl / 100 µl aqua bidest
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided

Only for research use

N:\WINWORD.DAT\IMMUNOCHEMIKALIEN\DATEN\BLAETTER\ANTIKOERPER\MISCELLANEOUS\AE1134.DOC