

Data Sheet

HUMAN PNEUMADIN

ANTIBODY, POLYCLONAL

Catalog no.:	AE1134.1 / AE1134.2
Immunogen:	Synthetic human Pneumadin (AGEPKLDAGV), poly Lysin conjugated
Synonyms:	PNM
Swiss-Prot No:	P22103
Gene Information:	Gene Name: N/A GeneID: N/A
Host:	Rabbit
Matrix:	Serum
Specificity:	Synthetic human Pneumadin
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	ELISA This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); -20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	12 January 2018

For research use only

Publishing research using AE1134? Please let us know so that we can cite your publication as a reference.

