



Data Sheet

HUMAN PNEUMADIN SYNTHETIC

Catalog no.:	AE1134AG.1 / AE1134AG.2
Name:	Human Pneumadin
Synonyms:	PNM
Sequence:	AGEPKLDAGV
MW:	0.96 kDa
Swiss-Prot No:	P22103
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	24 April 2018

For research use only

Publishing research using AE1134AG? Please let us know so that we can cite your publication as a reference.

