

## *DATA SHEET*

# ANTI HUMAN SECRETIN (aa 14-27) ANTIBODY, POLYCLONAL

<b>Catalog no.:</b>	AE 1135.1 / AE 1135.2
<b>Immunogen:</b>	Synthetic human Secretin (aa 14-27) poly Lysin conjugated (REGARLQRLLOGLV)
<b>Reference no.:</b>	LocusID 6343
<b>Host:</b>	Rabbit
<b>Matrix:</b>	Serum
<b>Specificity:</b>	Synthetic human Secretin (aa 14-27)
<b>Contents:</b>	20 µl / 100 µl (lyophilized) resuspend in 20 µl / 100 µl aqua bidest
<b>Known applications:</b>	ELISA (1: 5.00)

This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.

<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
------------------	--

### Only for research use

N:\WINWORD.DAT\IMMUNOCHEMIKALIEN\DATEN\BLAETTER\ANTIKOERPER\MISCELLANEOUS\AE1135.DOC