



Data Sheet

HUMAN SECRETIN (aa 14-27)

SYNTHETIC

Catalog no.:	AE1135AG.1 / AE1135AG.2
Name:	Human Secretin (aa 14-27)
Sequence:	REGARLQRLQLGLV
MW:	1.61 kDa
Swiss-Prot No:	P09683
Gene Information:	Gene Name: SCT Gene ID: 6343
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	20 April 2018

For research use only

Publishing research using AE1135AG? Please let us know so that we can cite your publication as a reference.

