

DATA SHEET

HUMAN SECRETIN (aa 14-27)

Catalog no.:	AE 1135AG.1 / AE 1135AG.2
Immunogen:	Human Secretin (aa 14-27) (REGARLQRLLOGLV)
Reference no.:	LocusID 6343
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in 0,1% acetic acid for required concentration
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided

Only for research use

N:\WINWORD\DAT\DATEN\BLETTER\ANTIKOERPER\MISCELLANEOUS\AE 1037.DOC