



Data Sheet

HUMAN PRO-ATRIAL NATRIURETIC PEPTIDE

[PRO-ANP] (aa 12-25)

ANTIBODY, MONOCLONAL

Catalog no.:	AE6201.1 / AE6201.2
Immunogen:	Synthetic human pro-ANP (aa 12-25) poly Lysin conjugated (LMDFKNLLDHLEEK(C))
Synonyms:	Cardiodilatin-related peptide
Swiss-Prot No:	P01160
Gene Information:	Gene Name: NPPA GeneID: 4878
Host:	Mouse Balb C
Clone no.:	A136-10 L13/10/H8
Isotype:	IgG ₁
Matrix:	Cell culture supernatant, Prot. G purified; 50 mM TRIS HCl pH 7,4
Specificity:	Synthetic human pro-ANP (aa 12-25)
Contents:	10 µg / 100 µg (lyophilized) Resuspend in 10 µl / 100 µl aqua bidest
Known applications:	ELISA <small>This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.</small>
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided

For research use only

