

Data Sheet

HUMAN PRO-ATRIAL NATRIURETIC PEPTIDE

[PRO-ANP] (aa 40-55)

ANTIBODY, MONOCLONAL

Catalog no.:	AE 6202.1 / AE 6202.2
Immunogen:	Synthetic human pro-ANP (aa 40-55) poly Lysine conjugated ((C)EPNEEAGAALSPLPEV)
Synonyms:	Atrial Natriuretic Factor, ANF
Swiss-Prot No:	P01160
Gene Information:	Gene Name: NPPA GeneID: 4878
Host:	Mouse Balb/c
Clone no.:	A137-M14-19-F9
Isotype:	IgG ₁
Matrix:	Cell culture supernatant, Prot. G purified; 50 mM TRIS HCl pH 7,4
Specificity:	Synthetic human pro-ANP (aa 40-55)
Contents:	10 µg / 100 µg (lyophilized) resuspend in 10 µl / 100 µl aqua bidest
Known applications:	ELISA <small>This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.</small>
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided

For research use only

