



Data Sheet

HUMAN pro-ATRIAL NATRIURETIC PEPTIDE [ANP] (aa 80-93)

SYNTHETIC

Catalog no.:	AE6202AG.1 / AE6202AG.2
Name:	Human pro-ANP (aa 80-93) (<u>C</u> SDRSALLKSKLRAL) Underlined amino acids differ from the original sequence; they were added or substituted for technical reasons.
Synonyms:	Atrial natriuretic factor (ANF)
Swiss-Prot No:	P01160
Gene Information:	Gene Name: NPPA, ANP, PND GeneID: 4878
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg/ 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided

For research use only



Immundiagnostik AG

Stubenwald-Allee 8a · 64625 Bensheim · Germany

Phone: +49 (0) 62 51/70 19 00 · Fax: +49 (0) 62 51/84 94 30 · info@immundiagnostik.com · www.immundiagnostik.com