



Data Sheet

HUMAN pro-ATRIAL NATRIURETIC PEPTIDE [ANP] (aa 80-93)

SYNTHETIC

Catalog no.:	AE6203AG.1 / AE6203AG.2
Name:	Human pro-ANP (aa 80-93) (SDRSALLKSKLRAL)
Synonyms:	Atrial natriuretic factor (ANF)
MW:	1557.7 Da
Swiss-Prot No:	P01160
Gene Information:	Gene Name: NPPA, ANP, PND GenelD: 4878
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg/ 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided

For research use only

