

DATA SHEET

HUMAN VASOACTIVE INTESTINAL PEPTIDE- RECEPTOR 2 [VIP-R2] (aa 3-14)

Catalog no.:	AE6302AG.1 / AE6302AG.2
Immunogen:	Human VIP-R2 (aa 3-14) (RFHLEIQEEETK)
Reference no.:	LocusID 7434
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in 0,1% acetic acid for required concentration
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided Repeated thawing and freezing must be avoided

Only for research use

N:\WINWORD\DAT\DATEN\LAETTER\ANTIKOERPER\CALCIUM AND BONE METABOLISM\A 7881.2.DOC