

DATA SHEET

HUMAN PARATHYROID HORMON [PTH]

(aa 1-38)

Catalog no.:	AK1105AG.1 / AK1105AG.2
Immunogen:	Human PTH (aa 1-38) (SVSEIQLMHNLGKHLNSMERVEWLRKKLODVHNFVALG)
Reference no.:	Locus ID 5741
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in 0,1% acetic acid for required concentration
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided

Only for research use

N:\WINWORD\DAT\DATENBLAETTER\ANTIKOERPER\CALCIUM AND BONE METABOLISM\IA 4223.2.DOC