



Data Sheet

HUMAN MYOSTATIN

RECOMBINANT

Catalog no.:	AK3003AG.1
Name:	Human Myostatin, full length
Synonyms:	Growth/differentiation factor 8 (GDF-8)
MW:	42.75 kDa
Swiss-Prot No:	014793
Gene Information:	Gene Name: MSTN, GDF8 GenelD: 2660
Source:	Produced serum-free in a human cell line
Purification:	Affinity chromatography
Storage buffer:	PBS pH 7.4
Purity:	> 95% (SDS-PAGE)
Contents:	10 µg (6.8 mg/ml)
Known applications:	ELISA, Western Blot This protein has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures.
Store at:	- 20 °C Repeated thawing and freezing must be avoided

For research use only



Immundiagnostik AG

Stubenwald-Allee 8a · 64625 Bensheim · Germany

Phone: +49 (0) 62 51/70 19 00 · Fax: +49 (0) 62 51/84 94 30 · info@immundiagnostik.com · www.immundiagnostik.com