

## Data Sheet

### HUMAN PRO-MYOSTATIN (aa 79-92)

#### ANTIBODY, POLYCLONAL

<b>Catalog no.:</b>	AK 3004.1 / AK 3004.2
<b>Immunogen:</b>	Synthetic human Pro-Myostatin (aa 79-92) KLH conjugated (SDGSLEDDDYHATT)
<b>Synonyms:</b>	Growth/differentiation factor 8 (GDF-8)
<b>Swiss-Prot No:</b>	O14793
<b>Gene Information:</b>	Gene name: MSTN, GDF8 GeneID: 2660
<b>Host:</b>	Rabbit
<b>Matrix:</b>	Serum
<b>Specificity:</b>	Human Pro-Myostatin (aa 79-92)
<b>Contents:</b>	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
<b>Known applications:</b>	ELISA (1:500)  This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided

#### For research use only

