



## DATA SHEET

# HUMAN PRO-MYOSTATIN (aa 79-92)

<b>Catalog no.:</b>	AK 3004AG.1 / AK 3004AG.2
<b>Peptide:</b>	Human Pro-Myostatin (aa 79-92) (SDGSLEDDDYHATT))
<b>Reference no.:</b>	LocusID 2660
<b>Source:</b>	Synthetic
<b>Purity:</b>	> 95% (HPLC)
<b>Contents:</b>	100 µg / 1 mg (lyophilized) Resuspend in 0,1% acetic acid for required concentration
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided

**Only for research use**

N:\WINWORD\DAT\DATENBLAETTER\ANTIKOERPER\CALCIUM AND BONE METABOLISM\A 4233.2.DOC