

## DATA SHEET

# ANTI HUMAN MYOSTATIN (aa 83-98) ANTIBODY, POLYCLONAL

<b>Catalog no.:</b>	AK 3005.1 / AK 3005.2
<b>Immunogen:</b>	Synthetic human Myostatin (aa 83-98) KLH conjugated (NMLYFNGKEQIIYGKI)
<b>Reference no.:</b>	LocusID 2660
<b>Host:</b>	Rabbit
<b>Matrix:</b>	Serum
<b>Specificity:</b>	Human Myostatin (aa 83-98)
<b>Contents:</b>	20 µl / 100 µl (lyophilized) resuspend in 20 µl / 100 µl aqua bidest
<b>Known applications:</b>	ELISA (1:500)
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided

**Only for research use**

N:\WINWORD\DAT\DATENBLAETTER\ANTIKOERPER\CALCIUM AND BONE METABOLISM\AK 3005.DOC