

DATA SHEET

HUMAN MYOSTATIN (aa 83-98)

Catalog no.:	AK 3005AG.1 / AK 3005AG.2
Peptide:	Human Myostatin (aa 83-98) (NMLYFNGKEQIIYGKI)
Reference no.:	LocusID 2660
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in 0,1% acetic acid for required concentration
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided

Only for research use

N:\WINWORD\DAT\DATENBLAETTER\ANTIKOERPER\CALCIUM AND BONE METABOLISM\A_4223_2.DOC