



Data Sheet

HUMAN MYOSTATIN (aa 53-66)

SYNTHETIC

Catalog no.: AK 3006AG.1 / AK 3006AG.2
Name: Human Myostatin (aa 53-66) (QKYPHTHLVHQANP)
Synonyms: Growth/differentiation factor 8 (GDF-8)
Swiss-Prot No: O14793
Gene Information: Gene name: MSTN, GDF8
GenelD: 2660
Source: Synthetic
Purity: > 95% (HPLC)
Contents: 100 µg / 1 mg (lyophilized)
Resuspend in 0.1% acetic acid for required concentration
Store at: 2-8 °C (lyophilized); - 20 °C (dissolved)
Repeated thawing and freezing must be avoided

For research use only

