



Data Sheet

HUMAN TISSUE TRANSGLUTAMINASE [TGc]

ANTIBODY, POLYCLONAL

Catalog no.:	AK 3010.1 / AK 3010.2
Immunogen:	Recombinant human TissueTransglutaminase [TGc]
Synonyms:	Protein-glutamine gamma-glutamyltransferase 2, EC=2.3.2.13, Transglutaminase C, Transglutaminase H, Transglutaminase-2
Swiss-Prot No:	P21980
Gene Information:	Gene Name: TGM2 GeneID: 7052
Host:	Rabbit
Matrix:	Serum
Specificity:	Human TissueTransglutaminase [TGc]
Contents:	20 µl / 100µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	ELISA (> 1:4000)
Store at:	2-8 °C (lyophilized); -20 °C (dissolved) Repeated thawing and freezing must be avoided

For research use only

