



DATA SHEET

HUMAN THYROGLOBULIN

Catalog no.:	AK3122AG.1 / AK3122AG.2
Protein:	Human Thyroglobulin
Reference no.:	LocusID 7038
Purity:	> 95% (SDS-PAGE)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in aqua bidest for required concentration
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided

Only for research use

N:\WINWORD\DAT\DATEN\BLETTER\PROTEINE\A 1005AG-proBNP1-76rec.DOC