



Data Sheet

MOUSE NIDOGEN 1

RECOMBINANT

Catalog no.:	AP1003AG.1 / AP1003AG.2
Name:	Mouse Nidogen 1
Synonyms:	Entactin
Swiss-Prot No:	P10493
Gene Information:	Gene Name: Nid1 GeneID: 18073
Purity:	> 95% (SDS-PAGE)
Matrix:	50 mM HEPES, 150 mM NaCl, pH 7.5
Contents:	10 µg / 100 µg (lyophilized) Resuspend in aqua bidest. for required concentration
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided

For research use only



Immundiagnostik AG

Stubenwald-Allee 8a · 64625 Bensheim · Germany

Phone: +49 (0) 62 51/70 19 00 · Fax: +49 (0) 62 51/84 94 30 · info@immundiagnostik.com · www.immundiagnostik.com