

Data Sheet

HUMAN CHROMOGRANIN A [CgA] (aa 284-301)

ANTIBODY, POLYCLONAL

| | |
|----------------------------|---|
| Catalog no.: | AST1007.1 / AST1007.2 |
| Immunogen: | Synthetic human CgA (aa 284-301) KLH-conjugated (<u>Y</u> QK <u>EEEEEM</u> AVVPQGLFRG) <small>Underlined amino acids were substituted or added for technical reasons.</small> |
| Synonyms: | Pituitary secretory protein I, SP-I |
| Swiss-Prot No: | P10645 |
| Gene Information: | Gene Name: CHGA GenelD: 1113 |
| Host: | Rabbit |
| Matrix: | Serum |
| Specificity: | Human CgA (aa 284-301) There was no cross reactivity detected with human CgA peptides (aa 1-17; 17-38; 63-76; 100-113; 116-130; 176-195; 238-247; 324-337; 361-372; 375-384; 411-424). |
| Contents: | 20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest. |
| Known applications: | RIA ^{1, 5, 6} , Western Blot, immunohistochemistry (paraffin sections, 1:80-1000) ^{1, 2, 3, 4, 6} <small>This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.</small> |
| Store at: | 2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided |
| References: | <ol style="list-style-type: none"> Stridsberg M, Oberg K, Li Q, Engstrom U & Lundqvist G (1995) Measurements of chromogranin A, chromogranin B (secretogranin I), chromogranin C (secretogranin II) and pancreastatin in plasma and urine from patients with carcinoid tumours and endocrine pancreatic tumours. <i>J. Endocrinol.</i> 144: 49–59. Portela-Gomes GM & Stridsberg M (2001) Selective processing of chromogranin A in the different islet cells in human pancreas. <i>J. Histochem. Cytochem.</i> 49: 483–490. Portela-Gomes GM & Stridsberg M (2002) Chromogranin A in the human gastrointestinal tract: an immunocytochemical study with region-specific antibodies. <i>J. Histochem. Cytochem.</i> 50: 1487–1492. Jin L, Zhang S, Bayliss J, Scheithauer B, Qian X, Kobayashi I, Stridsberg M & Lloyd R V (2003) Chromogranin a processing in human pituitary adenomas and carcinomas: analysis with region-specific antibodies. <i>Endocr. Pathol.</i> 14: 37–48. |



Antibodies

immundiagnostik

5. Stridsberg M, Eriksson B, Öberg K & Janson ET (2004) A panel of 11 region-specific radioimmunoassays for measurements of human chromogranin A. *Regul. Pept.* **117**: 219–227.

6. Tartaglia A, Portela-Gomes GM, Öberg K, Vezzadini P, Foschini MP & Stridsberg M (2006) Chromogranin A in gastric neuroendocrine tumours: An immunohistochemical and biochemical study with region-specific antibodies. *Virchows Arch.* **448**: 399–406.

Last updated on: 05 January 2016

For research use only

Publishing research using AST1007? Please let us know so that we can cite your publication as a reference.



Immundiagnostik AG

Stubenwald-Allee 8a · 64625 Bensheim · Germany

Phone: +49 (0) 62 51/70 19 00 · Fax: +49 (0) 62 51/84 94 30 · info@immundiagnostik.com · www.immundiagnostik.com