

Data Sheet

HUMAN CHROMOGRANIN B [CgB] (aa 480-493)

ANTIBODY, POLYCLONAL

Catalog no.:	AST1020.1 / AST1020.2
Immunogen:	Synthetic human CgB (aa 480-493) KLH-conjugated (<u>C</u> RFQDKQYSSHTAE) Underlined amino acids were substituted or added for technical reasons.
Synonyms:	Secretogranin-1, Chromogranin-B, CgB, Secretogranin I, Sgl
Swiss-Prot No:	P05060
Gene Information:	Gene Name: CHGB GenelD: 1114
Host:	Rabbit
Matrix:	Serum
Specificity:	Human CgB (aa 480-493)
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	RIA, Western Blot, immunohistochemistry (paraffin sections) This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	01 October 2015

For research use only

Publishing research using AST1020? Please let us know so that we can cite your publication as a reference.

