

DATA SHEET

ANTI HUMAN OXIDIZED FIBRINOGEN (INSOLUBLE FRACTION) ANTIBODY, POLYCLONAL

Catalog no.:	AT 1006.1 / AT 1006.2
Immunogen:	Human oxidized Fibrinogen (insoluble fraction)
Host:	Rabbit
Matrix:	Serum
Specificity:	Several oxidative modifications of Fibrinogen and Lipoproteins
Contents:	20 µl / 100 µl (lyophilized) resuspend in 20 µl / 100 µl aqua bidest

Known applications: ELISA (1: 500)

This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.

Store at: 2-8 °C (lyophilized); - 20 °C (dissolved)
Repeated thawing and freezing must be avoided

Only for research use

N:\WINWORD.DAT\DATEN\BLAETTER\ANTIKOERPER\OXIDATIVE STRESS\AT 1006.2.DOC