

## Data Sheet

# PORPHYROMONAS GINGIVALIS

## ANTIBODY, MONOCLONAL

<b>Catalog no.:</b>	AY1005.1 / AY1005.2
<b>Immunogen:</b>	Sonicated bacterial extract
<b>Host:</b>	Mouse Balb/c
<b>Clone no.:</b>	E3 P3 G4
<b>Isotype:</b>	IgG2b kappa
<b>Matrix:</b>	10 mM PBS pH 7.2, 150 mM NaCl, 0.099% NaN <sub>3</sub>
<b>Specificity:</b>	Porphyromonas gingivalis
<b>Contents:</b>	10 µg / 100 µg (liquid)
<b>Known applications:</b>	ELISA

This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.

**Store at:** 2-8°C (short term storage); -20°C (long term storage)

Repeated thawing and freezing must be avoided

**Last updated on:** 27 January 2016

### For research use only

**Publishing research using AY1005? Please let us know so that we can cite your publication as a reference.**

