

Data Sheet

PORPHYROMONAS GINGIVALIS

ANTIBODY, MONOCLONAL

Catalog no.:	AY1007.1 / AY1007.2
Immunogen:	Sonicated bacterial extract
Host:	Mouse Balb/c
Clone no.:	E5 P1 B8
Isotype:	IgG1 kappa
Matrix:	10 mM PBS pH 7.2, 150 mM NaCl, 0.099% NaN ₃
Specificity:	Porphyromonas gingivalis
Contents:	10 µg / 100 µg (liquid)
Known applications:	ELISA

This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.

Store at: 2-8°C (short term storage); -20°C (long term storage)
Repeated thawing and freezing must be avoided

Last updated on: 27 January 2016

For research use only

Publishing research using AY1007? Please let us know so that we can cite your publication as a reference.

