

Data Sheet

HUMAN CHOLINE KINASE ALPHA-1

ANTIBODY, MONOCLONAL

| | |
|----------------------------|---|
| Catalog no.: | AY1009.1 / AY1009.2 |
| Immunogen: | Recombinant human Choline Kinase alpha-1 |
| Synonyms: | Choline kinase alpha (CK), EC: 2.7.1.32, CHETK-alpha, Ethanolamine kinase (EK), EC: 2.7.1.82 |
| Swiss-Prot No: | P35790 |
| Gene Information: | Gene Name: CHKA, CHK, CKI GenelD: 1119 |
| Host: | Mouse Balb/c |
| Clone no.: | H-001 |
| Isotype: | IgM kappa |
| Matrix: | 10 mM PBS pH 7.2, 150 mM NaCl, 0.099% NaN ₃ |
| Specificity: | Human Choline Kinase alpha-1 |
| Contents: | 10 µg / 100 µg (liquid) |
| Known applications: | ELISA This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls. |
| Store at: | 2-8°C (short term storage); -20°C (long term storage) Repeated thawing and freezing must be avoided |
| Last updated on: | 29 January 2016 |

For research use only

Publishing research using AY1009? Please let us know so that we can cite your publication as a reference.

