

Data Sheet

HUMAN CHOLINE KINASE ALPHA-1

ANTIBODY, MONOCLONAL

Catalog no.:	AY1012.1 / AY1012.2
Immunogen:	Recombinant human Choline Kinase alpha-1
Synonyms:	Choline kinase alpha (CK), EC: 2.7.1.32, CHETK-alpha, Ethanolamine kinase (EK), EC: 2.7.1.82
Swiss-Prot No:	P35790
Gene Information:	Gene Name: CHKA, CHK, CKI GenelD: 1119
Host:	Mouse Balb/c
Clone no.:	H-006
Isotype:	IgM kappa
Matrix:	10 mM PBS pH 7.2, 150 mM NaCl, 0.099% NaN ₃
Specificity:	Human Choline Kinase alpha-1
Contents:	10 µg / 100 µg (liquid)
Known applications:	ELISA This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
Store at:	2-8°C (short term storage); -20°C (long term storage) Repeated thawing and freezing must be avoided
Last updated on:	29 January 2016

For research use only

Publishing research using AY1012? Please let us know so that we can cite your publication as a reference.

