

DATA SHEET

NORMAL MOUSE IGG

Catalog no.:	NG 4005.1/NG 4005.2
Host:	Mouse
Matrix:	Serum
Contents:	1 mg/5 mg mouse IgG fraction [NH ₄ (SO ₄) ₂ precipitation]
Store at:	2-8 °C (lyo.); - 20 °C (dissolved) Repeated thawing and freezing must be avoided

Only for research use

N:\WINWORD.DAT\IMMUNOCHEMIKALIEN\DATENBLAETTER\ANTIKOERPER\ANIMAL SERA\NG 4005.DOC