



Data Sheet

ANGIOTENSIN II

ANTIBODY, POLYCLONAL

Catalog no.:	A 1035.1 / A 1035.2
Immunogen:	Angiotensin II KLH-conjugated
Synonyms:	Angiotensin 2, Ang II
Swiss-Prot No:	P01019 (precursor Angiotensinogen is cleaved into Angiotensin-1, Angiotensin-2, and Angiotensin-3)
Gene Information:	Gene Name: AGT, SERPINA8 (Angiotensinogen) GeneID: 183 (Angiotensinogen)
Host:	Rabbit
Matrix:	Serum
Specificity:	Angiotensin II (100 % to human, mouse, rat, horse, dog)
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	ELISA This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	16 March 2022

For research use only

Publishing research using A 1035? Please let us know so that we can cite your publication as a reference.

